SCIENTIFIC REPORTS

Published online: 18 April 2019

OPEN Author Correction: Nucleus-specific X-ray stain for 3D virtual histology

Mark Müller¹, Melanie A. Kimm², Simone Ferstl¹, Sebastian Allner¹, Klaus Achterhold ¹ Julia Herzen¹, Franz Pfeiffer^{1,2} & Madleen Busse¹

Correction to: Scientific Reports https://doi.org/10.1038/s41598-018-36067-y, published online 14 December 2018

This Article contains errors.

In Figure 1A, the final step of the Hematein staining protocol is incorrectly labelled. The correct Figure 1 appears below.

Additionally, this Article contains a typographical error in the Results section, under the subheading 'MicroCT investigations of hematein stained mouse liver tissue' where,

"The tab water wash after staining is crucial for the specific cell nuclei staining result."

should read:

"The tap water wash after staining is crucial for the specific cell nuclei staining result."

This Article also contains a typographical error in the Methods section, under the subheading 'Hematein Staining Protocol' where,

"The soft-tissue sample was washed with 6 ml tab water (washing solution was changed every hour for the first 5 hours) and kept overnight in the washing solution."

should read:

"The soft-tissue sample was washed with 6 ml tap water (washing solution was changed every hour for the first 5 hours) and kept overnight in the washing solution."

Lastly, the Article contains an error in Table 1 where the incorrect % Hepatocytes value is given for Sub-cube 8. The correct value should read 27.0.

¹Department of Physics and Munich School of BioEngineering, Technical University of Munich, 85748, Garching, Germany. ²Department of Diagnostic and Interventional Radiology, Klinikum rechts der Isar, Technical University of Munich, 81675, Munich, Germany. Mark Müller and Melanie A. Kimm contributed equally. Correspondence and requests for materials should be addressed to M.M. (email: mark_mueller@ph.tum.de) or M.B. (email: madleen. busse@tum.de)




Figure 1. Staining protocol and interaction of the hematein-based X-ray stain with soft tissue. (**A**) The developed hematein-based staining procedure shows the individual steps involved including incubation and staining times. Staining step 1 was conducted using lead(II) acetate trihydrate as the heavy metal source. The lead(II) acetate trihydrate was dissolved in distilled water (c = 666 mM) and is referred to as working solution (A) (WS (A)). The staining step 2 involved a hematein solution in absolute ethanol (WS (B), 10% (w/v); c = 333 mM), which was derived from hematoxylin and was added to WS (A). More details concerning the entire staining protocol are described in the Materials and Methods section. (**B**) The positively charged hematein lead(II) complex (purple), which is built in situ in the soft-tissue sample, is interacting with the negatively charged phosphate backbone of the DNA (orange) present in the nucleus of the cell. The selective interaction of the hematein lead(II) complex with the DNA is achieved by acidification of the soft tissue during fixation or before staining and allows for a higher accumulation of the hematein lead(II) complex within the cell nucleus.

Open Access This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons license, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons license and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this license, visit http://creativecommons.org/licenses/by/4.0/.

© The Author(s) 2019