



Stable Isotope Dilution Analysis of the Major Prenylated Flavonoids Found in Beer, Hop Tea, and Hops

Lance Buckett¹, Simone Schinko¹, Corinna Urmann^{2,3}, Herbert Riepl^{2,3} and Michael Rychlik^{1*}

¹ Chair of Analytical Food Chemistry, Technical University of Munich, Freising, Germany, ² Organic-Analytical Chemistry, Weihenstephan-Triesdorf University of Applied Sciences, Straubing, Germany, ³ Campus Straubing for Biotechnology and Sustainability, Technical University of Munich, Straubing, Germany

Prenylated flavonoids from hops (Humulus lupulus) have become of interest in recent years due to a range of bioactivities. The potential health benefits of prenylated flavonoids include anti-cancerous activities and treatment of the metabolic syndrome among others. Since prenylated flavonoids from hops have shown pharmaceutical potential in clinical trials, robust analytical methods to determine their concentrations in food, supplements, and beverages are required. One such, the gold standard of analytical methods, is stable isotope dilution analysis due to its ability to compensate matrix effects and losses during sample work-up. As no commercial standards were available, the synthesis of seven different prenylated flavonoid isotopes utilizing various strategies (microwave assistance, acid base catalyst in the presence of deuterated substance and lastly, the use of Strykers catalyst) is described. The produced prenylated flavonoid isotopes were then applied in the first stable isotope dilution analysis method that quantified six natural prenylated flavonoids (Isoxanthohumol, Isoxanthohumol-C, 8-Prenylnaringenin, 6- Prenylnaringenin, Xanthohumol, and Xanthohumol-C) in beer, hop tea and hops to prove its applicability. The SIDA-LC-MS/MS method was validated resulting in LODs and LOQs for all analytes between 0.04 and 3.2 µg/L. Moreover, due to the simple clean-up the developed method allows the prospect for measuring clinical samples in the future.

Keywords: beer (analysis method for), stable isotope dilution analysis (SIDA), xanthohumol, xanthohumol C, isoxanthohumol, isoxanthohumol C, 8-prenylanaringeninx, 6-prenylnaringenin

INTRODUCTION

Prenylated flavonoids (PF) derived from hops have undergone an increased interest in recent years regarding their potential for improving health and uses in treating various disorders. Xanthohumol (XN), the most prevalent PF in hops shows an effect on metabolic syndrome (perhaps via modulation of glucagon -like-peptide-1), prevention of DNA damage, and has anti-cancer properties (1–4) Interestingly, many of XN plant metabolites and degradation products (**Figure 1**) possess bioactivities (4, 5). 8-Prenylnaringenin (8-PN) is the most potent phytoestrogen known so far (6). In addition 6-Prenylnaringenin (6-PN), the isomer of 8-PN, is slightly estrogenic, but also shows an anti-proliferative effect on different cancer cell lines (7). Another PF, Xanthohumol-C (XN-C), initiates differentiation in neuronal precursor cells and is neuroprotective (8). Accordingly, it seems to be a candidate for treating neurodegenerative diseases (9). The isomerized XN-C,

OPEN ACCESS

Edited by:

Sascha Rohn, Technische Universität Berlin, Germany

Reviewed by:

Michal Halagarda, Kraków University of Economics, Poland Giuseppe Annunziata, University of Naples Federico II, Italy

> *Correspondence: Michael Rychlik michael.rychlik@tum.de

Specialty section:

This article was submitted to Food Chemistry, a section of the journal Frontiers in Nutrition

Received: 21 October 2020 Accepted: 26 November 2020 Published: 15 December 2020

Citation:

Buckett L, Schinko S, Urmann C, Riepl H and Rychlik M (2020) Stable Isotope Dilution Analysis of the Major Prenylated Flavonoids Found in Beer, Hop Tea, and Hops. Front. Nutr. 7:619921. doi: 10.3389/fnut.2020.619921

1



isoxanthohumol-C (IXN-C) shows only a slight differentiation inducing effect and accordingly, that makes the chalcone structure all the more interesting (10). Clinical trials show that PFs have a pharmaceutical use as non-modified forms (flavones) or by slight modifications e.g., removal of the Michael system (i.e., cyclization of chalcones into flavanones, or hydrogenation of the Michael system) (11–13). Further studies investigating the *in-vivo* effects are required to fully understand the beneficial properties of PF concerning different diseases and, therefore, robust and sensitive analytical methods that target the bioactive hop compounds are needed, for example, stable isotope dilution analysis (SIDA) (1). To date, SIDA has not been applied for analysis of hop PFs.

Common analytical methods applied on PF found in hops are HPLC, immunoassays, square-wave adsorptive-stripping voltammetry, high performance thin layer chromatography, MALDI, and LC-MS/MS using two different ionization sources (APCI and ESI) (1, 14-16). Some of the caveats of the current methods are: (i) they require expensive development of antibodies, (ii) are not specific enough, (iii) rely on expensive SPE clean-up, (iv) have no internal standard, or (v) use a single internal standard which might have unseen errors due to matrix effects in MS measurements (1, 17). All the disadvantages of the methods above are overcome by SIDA, the current analytical gold standard in mass spectrometry (18). The advantages of SIDA compared to LC-MS/MS using an internal standard are that it can compensate the matrix effect and work up errors as only the ratio of the isotopolog is measured. The isotope labeling of the IS, guarantees accuracy also in sample workup due to having nearly the same physical and chemical properties as the analyte i.e., losses during workup are retained as the ratio of A/IS, which is commonly accepted to be constant (18). Combined with a simple clean up, SIDA is a rapid analysis method, which is suitable for large scale intervention trials. Therefore, we aim to develop the first SIDA method on hop PFs and test the appropriateness by quantifying PF content in beer, hops, and hop tea.

As the name suggests, SIDA requires the use of isotopelabeled substances. So far, XN was isotope-labeled via multi-step synthesis or plant feeding experiments. Initially, Berwanger et al. produced a ¹³C-XN by feeding hop flower cuttings with different labeled precursors and concluded that universally labeled ¹³C glucose was the best supplement with 9.41% of the product labeled (19). Additionally, the same research group produced ¹⁴C XN at 318 microCi/mmol (20). More recently, (21) carried out the successful total synthesis of XN enriched with 13 C (21, 22). The total synthesis of XN used seven steps producing a yield of 5.7%. Both of the above options are time consuming and expensive with low yields and no methods for other PFs other than XN are described. Furthermore, plant feeding experiments are generally poorly suited for SIDA, due to the presence of unlabeled isotopologues. In this study an approach to produce isotopically labeled PFs enriched with deuterium by utilizing various methods e.g., Strykers Cat, acid base catalyzation and microwave enrichment is described (23–27). The produced stable isotopes were then applied and validated using an LC-MS/MS-SIDA method to quantify the natural occurring PFs in a few selected samples of beer, hops, and hop tea.

MATERIALS AND METHODS

General Experimental Information

All reactions were carried under inert atmosphere of argon or nitrogen unless stated. All chemicals were of technical grade unless stated. TLC monitoring was carried out during all reactions. The TLC plates used for development were silica on alumina 60 F254 (Merck). For TLC development a solvent system of 60:40 EtOAc:Hexane was utilized (Merck, Deutschland) and the stain was prepared using CAM-reagent: Ammonium molybdate tetrahydrate (Sigma aldrich) (3 g), cerium ammonium sulfate dihydrate (Sigma aldrich) (1g), sulfuric acid 10% 100 mL (Sigma aldrich), stirred vigorously and filtered twice through filter paper (H 1/2 filter paper Schleicher & Schuell). TLC plates were dipped in the stain and heated using a heat gun until the TLC plates turned blue and dark spots formed. The analytical standards that underwent qNMR were dissolved in 600 µL MeOD-d4 (≥99.8% D, Sigma Aldrich, USA) with NMR tubes (177.8 \times 4.97 mm, Bruker, USA). Staff members at the Chair of Molecular Sensory Science at the Technical University Munich carried out the NMR and qNMR experiments with a Bruker Advance II 400 MHz and 500 MHz NMR spectrometer. Data evaluation was conducted on MestReNova (Mestrelab Research, Spain) software. All samples that underwent LC-MS/MS analysis were filtered using a 0.22 μm filter (Chromafil[®] Xtra PVDF-20/13, 0.20 μm, Ø 13 mm; Macherey-Nagel).

LC-MS/MS Analysis

Preparation of Standard Solutions

Each analytical and internal standard underwent quantitative NMR (qNMR) in methanol-d6 to determine the concentration before being made up to a final concentration of 1 mg/mL and placed in -80° C freezer. An exception is D3-8-PN and D3-XN-C the concentrations of which were calculated using HPLC-UV from calibration graphs produced by the standard compounds 8-PN and XN-C as the amounts produced in the reactions were too low for qNMR.

LC-MS/MS Parameters

The chromatographic separation of the analytes was performed on a Shimadzu Nexera X2 UHPLC system (Shimadzu, Kyoto, Japan). The LC-parameters used solvent (A) H2O (0.1% FA) with ACN (0.1%) FA, solvent (B) with the initial solvent settings started at 25% B and increased to 35% B at 4 min where it remained until being linearly increased from 35 to 95% at 12.5 min. At 12.5 min the solvent remained at 95% B for 3.5 min until returning over 1 min back to 25%. Overall, each run took 22.5 min. The MS instrument operated in the positive electrospray ionization (ESI) mode for all analytes. The column (Restek Raptor FluoroPhenyl 2.7 µm 100 x 2.1 mm) protected with a precolumn with penta fluoro propyl phenol (PFPP) Shimadzu, Velox EXP Guard 2.1 x 5 mm) was equilibrated for 5 min between each injection. Samples were injected at 1 μ L, the flow was set at 0.3 μ L/min and the column oven was set to 30°C. Using an additional valve after the column, the solvent flow was introduced into the mass spectrometer at 4 min after injection to prevent the instrument from matrix contamination. The LC was coupled to a triple quadrupole mass spectrometer (LCMS-8050, Shimadzu Corporation, Kyoto, Japan). Parameters for MS were as follows: heating gas at a flow 10 L/min, nebulizing gas flow 3 L/min, drying gas flow 10 L/min, heat block temperature 400°C, desolvation line temperature 250°C, interface temperature 300°C, interface voltage 4 kV, and collision-induced dissociation gas pressure 270 kPa. The mass spectrometer worked in the scheduled multiple reaction monitoring (MRM) mode for MS/MS measurements. Optimized voltages and collision energies, the retention time of each analyte are listed in Table 1. The LabSolutions software (Shimadzu, Kyoto, Japan) was utilized for data acquisition and data analysis.

Analytical HPLC

The analytical HPLC analysis utilized a Shimadzu HPLC-DAD (SPA-M20A) coupled with a liquid chromatograph system (LC-20 AD, Shimadzu, Kyoto, Japan). Samples were automatically injected with an autosampler (SIL-20A HAT, Shimadzu, Kyoto, Japan) and a degassing unit (DGU-20A3R) was used when solvent system was active. Analyte separation used a C18 (YMC-Pro pack, 3μ m, 12 n, 150×3 mm Machery-Nagel) with a precolumn (C18, 8×3 mm). The mobile phase consisted of 0.1% FA in HPLC grade water as solvent A and 0.1% FA in ACN as solvent B. The flow rate was set 0.3 ml/min with gradient starting at 50% B for 2.0 min, and then gradually increasing solvent B to 90% at 11.0 min. The mobile phase remained at 90% for a further 2.5 min before being returned to 50% B. The column was

equilibrated for 5 min between injections of 10 μ L. Chalcones were monitored at 370 nm and flavanones at 290 nm. 8-PN was used for the calibration of D3-8-PN and XN-C was used for D3-XN-C at 0.56, 1.1, 1.7, 2.2, 2.8, 3.3, 3.9, 4.4, and 5.0 mg/L in ACN (see **Supplementary Material**).

Sample Preparation Beer

Fifteen beers that represented a cross section of brewery products of Bavaria were purchased from local supermarkets. They were opened and degassed using sonication. A 1 mL aliquot of each beer was placed into a 4 mL vial and IS were spiked at various concentrations, deuterated isoxanthohumol (D3-IXN) at 6.99 μ M, deuterated isoxanthohumol-C (D3-IXN-C), and deuterated xanthohumol (D3-XN-C) at 0.14 μ M, deuterated 8-prenylnaringenin (D3-8PN) at 0.29 μ M, deuterated 6-prenylnaringenin (D3-6-PN) at 0.15 μ M and deuterated xanthohumol (D3-XN) at 0.14 μ M. 1 mL of EtOAc was added and the vials mixed using a vortex (Vortex Genie 2, Scientific Industries). Taking the supernatant, which was repeated 3 times, the combined organic phase of the samples was dried under N2 (set at 30°C) and resolved in 1 mL ACN.

Hops

Using a method derived from Stevens et al., a packet of pelletized hops (Columbus, 15.6% α -acid, USA HOPS) of a high alpha acid (15.6% AA) was purchased from the local brewing shop and ground frozen using a mortar and pestle (28). 0.1 g of the powder was taken and extracted using sonication and 100 mL of EtOAc for 15 min. Samples (0.1 mL) were spiked with IS, dried and re-dissolved in 1.0 mL of ACN.

Hop Tea (5%)

A brand of hop tea was purchased from the local supermarket and two different preparations were investigated. One tea bag was prepared following the instruction on the packet. 100 mL of boiling (distilled) water was added to one tea bag and brewed for 5 min. Once the tea was at room temperature the IS were spiked and the solution was dried using a N2 drier and a heating plate set at 30°C. Another sample of tea was crushed using a mortar and pestle and prepared the same as the hop pellets.

Method Validation and Statistical Analysis

The developed method and workup procedures were validated according to the procedure of Vogelgesang and Hädrich (29). All LC-MS/MS data was collected and integrated using lab solutions software (Shimadzu Europa GmbH) and data transferred into Microsoft excel where all calculations were performed using excel functions (Microsoft Corporation, USA).

Preparation of Calibration Graphs

Calibrators of all 6 analytes (Reference standards) were spiked against their internal standards (IS) in ACN at the following molar ratios [n(Ref)/n(IS)]: IXN, 14 points: 100:1, 80:1, 60:1, 40:1, 20:1, 10:1, 7.50:1, 1:1, 1:20, 1:40, 1:60, 1:80, 1:100 ($R^2 = 0.9982$). IXN-C 12 points: 100:1, 60:1, 25:1, 10:1, 7.5:1, 5:1,1:1, 1:20, 1:40, 1:60, 1:80, 1:100. ($R^2 = 1$). 8-PN, 14 points: 25:1, 15:1, 12.5:1, 9.7:1, 7.5:1, 5.0:1, 4.4:1, 3.9:1, 3.3:1, 2.8:1 1.7:1,

Analyte-D ₃	D ₃ -XN	D ₃ -IXN	D ₃ -6-PN	D ₃ -8-PN	D ₃ -IXN-C	D ₃ -XN-C
MRM-1	358.20–180.15	358.25–180.05	345.20–289.20	345.20–166.15	356.20–235.20	356.20–234.20
	(CE:–21.0)	(CE:–26.0)	(CE:–16.0)	(CE:–25.0)	(CE:–22.0)	(CE:–20.0)
MRM-2	358.20-302.20	358.25–302.15	345.20–167.10	345.20–167.15	356.20–234.70	356.20–81.10
	(CE:-12.0)	(CE:–17.0)	(CE:–27.0)	(CE:–24.0)	(CE:–23)	(CE:–16.0)
MRM-3	358.20–181.15	358.–301.20	345.20–166.15	345.20–288.20	356.20-82.10	356.20–111.15
	(CE:–23.0)	(CE:–15.0)	(CE:–28.0)	(CE:–16.0)	(CE:-44.0)	(CE:–33)
Analyte-Ref	XN	IXN	6PN	8-PN	IXN-C	XN-C
MRM-1	355.10–179.25	355.30–179.25	341.35–165.20	341.35–165.20	353.00–233.20	353.20–233.15
	(CE:–21.0)	(CE:–24.0)	(CE:–29.0)	(CE:–18.0)	(CE:–22.0)	(CE:–21.0)
MRM-2	355.10-299.30	355.30-299.30	341.35–285.30	341.35–285.25	353.00–191.20	353.20–191.10
	(CE:-12.0)	(CE:-15.0)	(CE:–15.0)	(CE:–15.0)	(CE:–32.0)	(CE:–31.0)
MRM-3	355.10-113.25	355.30-172.15	341.35–123.10	341.35–123.0	353.00-109.20	353.20–109.15
	(CE:-23.0)	(CE:-9.0)	(CE:–37.0)	(CE:–33.0)	(CE:-35.0)	(CE–33.0)
Deuterium Enrichment	m/z +3	m/z + 3	m/z +4	m/z +4	m/z +3	m/z +3

TABLE 1 | The MRM profiles of the isotopologues compared with the MRM profiles of the reference compounds being analyzed.

0.1:1, 0.05:1, 0.04:1 ($R^2 = 0.9995$). 6-PN 11 points: 20:1, 15:1, 12.5:1, 5.9:1, 1.98:1, 1.5:1, 1.3:1 0.65:1, 0.43:1, 0.2:1, 0.002:1. (R = 0.9999). XN, 11 points: 100:1, 80:1, 60:1, 40:1, 20:1, 1:1, 1:20, 1:40, 1:60, 1:80, 1:100 ($R^2 = 0.9996$). XN-C, 10 points: 80:1, 60:1, 40:1, 20:1, 1:1, 1:20, 1:40, 1:60, 1:80, 1:100 (R^2 0.9997). The calibration graphs were formed from the LC-MS/MS data collected and from integration of the area of each peak using lab solutions software (Shimadzu Europa GmbH). The area of each peak A(ref)/A(IS) was then plotted into Microsoft excel (Microsoft Corporation, USA) against the known n(ref)/n(IS). Mandel testing confirmed linearity of the graphs and the equation resulted in determining the concentration of the prenylated flavonoids found in the samples (30). See Supplementary Material for each plotted calibration graph. Each calibrant level was injected 3 times and are the average value reported.

Precession, Limits of Detection and Limits of Quantification

The LOD concentrations were initially estimated by DIN 32645 (31). A hop free beer, a kindly gift of the TUM Research Center Weihenstephan for Brewing and Food Quality., served as blank matrix. Activated charcoal 1 g was added to the hop free beer (50 mL) and the solution was left to stirring for 1 h before being filtered through filter paper. The hop free beer was extracted three times with EtOAc (50 mL) and then concentrated to dryness and dissolved in ACN. The blank matrix was then spiked with the corresponding amounts of analytes near the estimated LOD to apply the method by Vogelgesang and Hädrich (29): IXN (0.18, 0.70, 1.3, 1.8 µg/L, and IS at 1 µg/L), IXN-C (0.05, 0.35, 0.2, 0.5 μ g/L, and IS at 0.5 μ g/L), 8-PN and XN (0.25, 1.0, 1.8, 2.5 μ g/L, and IS at 2.0 µg/L), 6-PN (0.16, 0.64, 1.1, 1.6 µg/L, and IS at 1.0 μ g/L), XN-C (0.2 μ g/L, 0.8 μ g/L, 1.4 μ g/L, 2.0 μ g/L and IS at 0.5 μ g/L). For recoveries the matrix was spiked at 4 different levels (2 low and 2 high). IXN low level, 0.004, 0.009 mg/L (D3-IXN 0.0025 mg/L), and high 2.5, 4.0 mg/L (D₃-IXN at 2.5 mg/L). IXN-C low level, 0.004, 0.009 mg/L (D₃-IXN-C 0.0025 mg/L), high level 0.5, 0.05 mg/L (D₃-IXN-C 0.1 mg/L). 8-PN, low level 0.004, 0.009 mg/L (D₃-8-PN, 0.0025 mg/L). High level 0.5, 0.25 (D₃-8PN, 0.1 mg/L). 6-PN low level 0.004, 0.009 mg/L (D₃-6-PN 0.00125 mg/L). High level 0.5, 1.0 mg/L (D₃-6PN 0.25 mg/L). XN, low level 0.004, 0.009 mg/L (D3-XN 0.0025), and high level 0.5, 1.0 mg/L (D₃-XN 0.5 mg/L) XN-C low level, 0.004 mg/L and 0.009 mg/L (D₃-XN-C: 0.0025 mg/L), and high level 0.025, 0.05 (D₃-XN-C 0.05 mg/L). Each sample underwent the same extraction in section Beer, and molarity of each compound was calculated using linear regression using the calibrations graphs prepared above. The ratios of the theoretical spiked vs. the calculated values determined the recoveries, LOD and LOQ. Precision was checked by triplicate determination of a mixed sample being injected 9 times and the RSD of ratios of the A(ref)/A(IS) calculated. Intradaily precision was determined by measuring 3 samples twice during one day and the inter-daily precision was determined by measuring a sample 3 times weekly using double injections.

Calculation of the Prenylated Flavonoid Content

The concentration of each analyte in each sample was calculated by taking the integrated area ratio of the analyzed reference compound (unknown) compared to the peak area of the spiked internal standard. The ratio of A(ref)/A(IS) was then used with linear regression equation reported in each excel graph to calculate the concentration [i.e., n(unknown)/n(IS)] in the samples in triplicate +- SD. The total PF content was calculated from the sum of each analyte using the excel function. The concentration of the PF in samples are expressed as the mean +- SD.

Analytical HPLC

To characterize and quantify labeled XN-C and 8-PN analytical HPLC was necessary. HPLC using a LC-20 AD Shimadzu Prominence HPLC System (Shidmadzu Kyoto, Japan) coupled with a SPA-M20A prominence diode array detector. And a YMC-

Pack Pro C18, S-3 μ m, 12 n, 150 x 3 mml. D, precolumn: 8 x 3 mm (Machery-Nagel, Düren, Germany) column. The retention times and UV profiles were compared with the non-deuterated compound (see **Supplementary Material**). The non-deuterated forms of XN-C and 8-PN were constructed into calibration solutions of 5.0, 4.4, 3.9, 3.3, 2.8, 2.2, 1.7, 1.1, and 0.6 mg/L. Upon evaluation of the graphs equations of XN-C y = 0.00001x + 0.0724, R^2 = 0.998 and 8-PN y = 42810x + 1007, R^2 = 0.999 were used to quantify the concentration of the D3-XN-C and D3 8-PN.

RESULTS AND DISCUSSION

Synthetic Results

As no commercial isotopically labeled standard of any prenylated hop flavonoid were available all isotopologues were synthesized (Figure 2) before SIDA was carried out (for detailed synthesis instruction see Supplementary Material). First, the isotope labeled standard of IXN was synthesized using wet chemistry, cyclisation following (32), with all protic solvents replaced by deuterated equivalences (32). To produce isotopically labeled xanthohumol (D₃-XN), 6-prenylnaringenin (D₃-6PN), 8-prenylnaringenin (D₃-8PN) and isoxanthohumol (D₃-IXN) two synthetic methods were combined. One, the removal of aryl methyl ethers by LiCl in DMF and the other, enrichment of deuterium by D₂O and catalyst mixtures of Pt (30%)/C, Pd (30%)/C, Rh (10%)/Al. (26, 33). Using this approach, four isotopically labeled standards (D3-IXN, D3-6-PN, D3, 8-PN, and D₃-XN) were created in 1 pot (Supplementary Material) reaction. A recent method using the demethylation and microwave approach has greatly improved the yields for production of non-deuterated 6- and 8-PN (27). The synthesized D₃-XN and deuterated IXN (by wet chemistry approach) were then further reacted into D₃-XN-C and D₃-IXN-C using a known selective cyclisation of the prenyl-group (10).

To overwhelm the electronically activation of the Michael system of XN, Stryker catalyst, which selectively lead to hydrogenation of the double bond, was used to synthesize the isotope labeled standard dihydroxanthohumol (DXN) by using a deuterated form rather than conventional donation of protons (**Supplementary Material**) (23, 34).

Mass Spectrometry of Synthesized Isotopes

The enrichment of deuterium in the isotopologues were confirmed using MS (see **Table 1**) and NMR analysis (**Supplementary Material**). In MS analysis, PFs fragmentate mainly via retro-diels alder which leads in case of the natural chalcone XN to the A-ring fragment (m/z 179) and the B-ring fragment (m/z 120). The fragmentation pattern of the isotopolog shows an A-ring (m/z 180), which leads to the assumption, that one deuterium is exchanged on the A-Ring. Furthermore, it can be concluded from fragmentation pattern, that one deuterium is exchanged on the B-Ring and one on the Michael system (**Figure 3**). Importantly no labeling is in the prenyl side chain as this would not be visible in MRM fragmentation and the fragment ion would have the same m/z as the precursor ion. A

similar pattern can be seen in the analysis of D₃-IXN, D₃-6PN, and D₃-8PN (**Table 1**). In case of the DXN derivative the high resolution mass spectrometry showed the incorporation of only one deuterium on the former α , β -unsaturated double bond (**Supplementary Material**), leading to the assumption that keto-enol equilibration does not save the deuterium in α -position from being exchanged. The NMR analyses confirmed the results of the mass spectrometry and are shown in the **Supplementary Material**.

The SIDA method should allow a simple, rapid and robust analysis concerning prenylflavonoids (PFs) for quantification in beer, hop products and future in vivo experiments. For the latter the DXN derivative is suitable, since the Michael system suggests a reduced bioavailability. All substances have a natural distribution of isotopes (in particular carbon-12 and carbon-13) and for an accurate SIDA there needs to be a certain increase in the mass making the difference measurable by mass spectrometer. The fact that DXN only has one D label made it not viable to be used in the SIDA method. However, the isotope analogs of the main natural prenylflavonoids D₃-IXN, D₃-IXN-C, D₃-8-PN, D₃-6PN, D₃-XN, and D₃-XN-C were deemed viable as they all had an increase between m/z of +3 to +4 Daltons. Table 1 shows a summary of the results for collisions energies and the mass transitions that were selected during the MRM optimizations and were incorporated into the LC-MS/MS method (35).

Method Validation

Calibration for SIDA was performed by measuring mixtures of the reference standards and the internal standards at known concentrations. The response factor was obtained by linear regression plots after testing linearity according to (30). The calibration curves for each compound plotted are Area(ref)/Area(IS) over the n(ref)/n(nIS) and resulted from a triple injection (Supplementary Figures). The LOQ and LOD was calculated in Microsoft Excel (Microsoft Corporation, USA) according to the requirements of Vogelgesang and Hadrich, i.e., blank matrices were spiked with reference compounds and internal standards at known concentrations near an estimated LOD (29). The samples were worked up according to section Beer and analyzed in triplicate and concentrations of each analyte were compared to the theoretical values. The LOD was set at the 95% confidence level, where the maximum peak of a blank would be detected (29). The LOQ was set at four times the LOD. In a sense the LOQ and LOQ was a recovery test, but performed at an extreme. Therefore, additional recoveries, were tested at different concentrations of analytes and not above the LOQ. The recoveries were determined on a low (0.004-0.009 mg/L) and a high level (0.025-4.0 mg/L) resulting in 81.0-113.3% (see Table 2) for each specific analyte.

The LOD and LOQ (**Table 2**) are similar to previously reported values, although the levels for XN and IXN presented here are lower. Stevens et al., using a triple quadrupole LC-MS/MS system with selective reaction monitoring, reached a LOQ of 10 ng/mL, which is at least double of the values measured here (LOQ 0.3–3.74 ng/mL) (28). The improvement is probably due to the use of SIDA with the carrier effect using



internal standards eluting simultaneously with the analytes (28). However, using the more conventional method, i.e., comparing the signal to noise ratio, might provide a lower LOD than the method of Vogelgesang, and Hadrich, as it has less stringent criteria by not spiking blank matrixes and reporting the response (29). Using a single quadrupole MS PFs were analyzed in urine with LODs between 0.2 and 0.6 ng/mL (36). Additionally, in the method described here only 1 μ L of sample was injected into the MS and therefore, room for improvement might be possible as increasing the injection volume might increase the sensitivity. Perhaps for the analysis of clinical samples this would be necessary.

Quantification of Hop Prenylated Flavonoids in Samples Using SIDA

For SIDA a simple LLE method was developed. It enabled a faster work up than solid phase extraction (SPE). Furthermore, it has been reported that LLE is superior in regards to retaining the PF as SPE results in poor recoveries for some analytes as each PF is eluted at different retention factors (1). This was also observed by us using SPE in the development of the analytical method. Spiking LLE samples with labeled standards allowed minimal clean-up steps and less losses of analytes due to the ratio of each analyte against the internal standard remaining constant throughout the work-up process as the physical properties of the internal standard and the analyte are virtually identical. Therefore, regardless of the total losses of analyte the concentration can be calculated as long as the analyte and the internal standard ratio is detectable. The matrix effect during beer analysis was anticipated by utilizing SIDA due to the respective labeled internal standards being eluted at the same time. A structurally different internal standard has a different elution time than the analyte, which can lead to large errors as there is no compensation for the matrix suppression or enhancement that may occur. Additionally, here for the first time we report the quantitation of IXN-C and XN-C in beer, which are relatively unstudied prenylated flavonoids. Lastly, it is the first SIDA method that has been utilized on prenylated flavonoids found in hops.

LC-Parameters

The LC method was designed so that sugars and carbohydrates that may have been extracted during the LLE were mostly omitted by diverting the initial 6 min of LC to waste. Additionally, the elution of the compounds of interest was devised to the isocratic phase of the method enabling a more uniform peak shape. The PFPP LC column used has been reported by Sus et al. (37) to be very suitable to separate prenylated flavonoids due to the phenolic-phenolic interactions (37). Additionally, compared to the C₁₈ phase, no apparent retention time shifts were observed with the PFPP stationary phase (38). Thus, most of the analytes are eluted separately and IXN-C and 8-PN are almost base line separated (**Figure 4**) which is not a problem in MS quantification.

Quantification of Prenylated Flavonoids in Beer, Hop Tea, and Hops

In order to prove the applicability of the new SIDA we quantified the PFs in different beers as summarized in **Table 3**. As the selection of beer samples are neither representative nor comprehensive, our results may only indicate the principle distribution of PFs in the different beer types. The total concentrations of PFs reported by (28) were between 43 and 4000 μ g/L, and in European beers between 43 and 2680 μ g/L (28, 39). Here we report comparable results between 723 and 2456 μ g/L (see **Table 3**). A more recent analysis of beer polyphenols revealed a range of prenylated flavonoids of 0 and 9500 μ g/L and highest containing beers were stouts and IPAs with an average range between 1.900 and of 2190 μ g/L (40). The values of the stouts in general are comparable to



our values (2445 μ g/L), and only one sample (a black IPA with 9500 μ g/L PF) which was not analyzed in our study, contained vastly more PF. It is worth mentioning that the alcohol free beer has a significantly reduced prenylated flavonoids content (1350 μ g/L) compared to the IPA with alcohol (2400 μ g/L). Obviously, the process of removing alcohol from beer seems to remove total PF content and reduce the amounts of compounds in the chalcone form, which are considered as desirable.

The amount of XN-C and IXN-C is very low in all beer samples, but it is shown that IXN-C is higher than XN-C in all beers measured, except for the double IPA. This is not surprising, since most of the XN is also converted to IXN during the brewing process (41). Interestingly, the amount of XN-C in the coffee stout is similar to the amount of IXN-C, which could suggest a stabilization process. This could be due to Maillard reaction products, which are supposed to stabilize XN from cyclizing in dark beer (42, 43). The double IPA is a heavily dry hopped beer and that might explain the higher XN-C amounts compared to IXN-C as the hops are not boiled to the extent of other beers. The chromanlike flavonoids XN-C and IXN-C are relatively unstudied and might confirm that they may add to the hepatoprotective properties that are said to outweigh the ethanol effects from beer consumption (44). However, this would need further research as little is known regarding the bioactivities of XN-C and IXN-C right now.

Furthermore, the concentrations of PFs in hop tea was investigated. However, the concentrations are so minimal that it would be very surprising if a health effect or taste modulation occurs (**Table 4**) (45). Drinking a tea that is brewed would only have 6.5 μ g PF (based on 300 mL). Comparing this to the consumption of an alcoholic free beer which would contain 440 μ g is vastly different. One optimistic observation

Analyte	LOD µg/L	LOQ µg/L	Recovery% High 2-points $(n = 3)$	Recovery % Low 2-points $(n = 3)$	Intra-day Precision RSD% (<i>n</i> = 3)	Inter-day Precision RSD% ($n = 3$)	Injection Precision RSD% (n = 9)
IXN	0.07	0.30	89.4–105	95.0-104	3.3	6.0	0.51
IXN-C	0.21	0.84	89.9-102	93.6-113	4.1	4.5	4.2
8-PN	0.32	1.3	106-109	108–113	2.09	8.2	4.87
6-PN	0.94	3.74	89.7-104	91.0-101	1.52	3.12	3.42
XN	0.06	0.24	90.0-109	90.9–93.0	3.15	0.51	1.24
XN-C	0.19	0.78	82.8–93.8	81.0-91.6	2.43	0.45	2.77

TABLE 2 | The LOD, LOQ recovery and precision determined in the hop free beer matrix of each analyte quantified using the developed SIDA method.



is that IXN concentrations in brewed tea is not increased, which might imply that using hot distilled water does not cause cyclisation to the less desired IXN, contrary to boiling hops, a process essential in the manufacture of beer, where most XN is thermally cyclized into IXN (41). Perhaps a formulation of hop tea comprised from hop pellets or XN extract might be an alternative supplement as they contain very high concentrations of prenylated flavonoids of 11.200 mg/kg reported here with the main compound being XN (6640). Although the solubility of XN in hot water is not high (46, 47). However, the solubilization of other (e.g., IXN, 6PN, 8PN) PF at room temperature in water is higher, which might make tea formulations of hops a viable, but extremely bitter tasting dietary source (47).

CONCLUDING REMARKS

The synthesis of the deuterated labeled internal standards applicable for prenylated flavonoids analysis has been developed

using three different methods. Firstly using a microwave assisted demethylation while simultaneously enriching compounds with deuterium was envisaged and applied. Second, utilizing functional groups such as a Michael system and Strykers Cat incorporating deuterium and lastly, we have also shown that a method does not need a microwave reaction, which some laboratories could be limited by a lack of resources. By contrast, acid and base reactions are also suitable to introduce deuterium. Using the synthesized products we conducted a robust SIDA method for the analysis of six different PF found in beer and hops with comparable results (723-2456 µg/L) to other investigations. For a statistical differentiation of the beer types a broader and more comprehensive sampling would be necessary. Regarding the investigation of the PF content in hop tea it contains extremely low concentrations of PF, i.e., it is very unlikely the content is contributing to health or taste (45). The method here uses a fast clean-up of LLE of samples containing PF and compensates for losses during workup by the addition of isotopically labeled internal standards. The LOD and LOQ are excellent and precision is within the DIN 32645

TABLE 3	Concentration of prer	nvlated flavonoids found in	15 beer styles using the deve	eloped SIDA method: \pm (SD μ a/L	n = 3
IADLE 3	Concentration of pref	nyialed liavonoids lound in	To been styles using the deve	eloped SIDA method: \pm (SD μ g/L	., // =

	IXN μg/L	IXN-C µg/L	8-PN μg/L	6-PN μg/L	XN μg/L	XN-C μg/L	Total μg/L
1. Helles (pale lager)	1160 (36)	12.0 (0.73)	19.7 (2.29)	46.7 (0.4)	19.8 (1.7)	1.9 (0.11)	1260
2. Festbier (festival beer)	1110 (26)	25.6 (6.30)	25.3 (1.14)	51.8 (0.90)	19.2 (2.5)	2.3 (0.09)	1234
3. Lager	627 (2.5)	13.7 (0.22)	14.4 (0.58)	14.2 (0.1)	22.1 (1.4)	1.9 (0.10)	693
4. Pilsner beer	989 (4.8)	20.5 (0.74)	7.0 (0.28)	22.8 (0.2)	10.1 (0.7)	1.8 (0.07)	1051
5. Hopped wheat beer A	1260 (53)	28.7 (0.74)	42.0 (1.07)	278 (5.6)	425 (35.1)	8.3 (0.45)	2042
6. Hopped wheat beer B	1400 (93)	15.4 (1.07)	36.5 (1.20)	120 (6.7)	411 (18.2)	10.2 (0.69)	1993
7. Double IPA	1320 (14)	12.8 (1.08)	50.6 (2.24)	421 (9.7)	500 (22.5)	16.1 (1.19)	2331
8. IPA	1540 (280)	14.1 (0.24)	54.6 (2.59)	368 (15.9)	464 (49.4)	5.1 (0.39)	2446
9. IPA alcohol free	1150 (147)	21.7 (4.96)	31.5 (1.98)	104 (3.2)	21.7 (3.9)	2.9 (0.37)	1332
10. Dunkel (dark beer)	1170 (128)	14.9 (1.35)	101.0 (7.55)	355 (3.1)	728 (30.0)	5.6 (0.20)	2375
11. Coffee stout	1480 (51)	23.4 (0.47)	88.7 (4.12)	351 (11.2)	483 (4.8)	18.8 (0.79)	2445
12. Bock beer	1550 (39)	19.0 (0.69)	50.3 (2.49)	156 (3.0)	99.3 (4.3)	4.4 (0.50)	1879
13. Doppelbock beer A	915 (15)	17.3 (1.75)	47.7 (2.64)	203 (2.40)	155 (29.9)	7.6 (0.09)	1346
14. Doppelbock beer B	1460 (80.6)	25.4 (1.67)	35.5 (1.45)	63.6 (0.6)	27.5 (1.8)	2.1 (0.25)	1614
15. Wheat bock beer	579 (13.9)	15.5 (1.00)	22.5 (1.12)	91.2 (3.2)	59.5 (3.9)	3.2 (0.20)	771

TABLE 4 | The concentration of the prenylated flavonoids in hops and hop tea \pm (SD $\mu\text{g/L}).$

IXN	IXN-C	8-PN	6-PN	XN	XN-C	Total
4250	151	4.46	75.5	6640	76.9	11,200
(274)	(1.86)	(0.16)	(11.4)	(78.2)	(1.54)	
14.4 (2.44)	28.1 (0.52)	5.93 (0.52)	<lod< td=""><td>17.9 (0.49)</td><td>3.12 (0.13)</td><td>69.0</td></lod<>	17.9 (0.49)	3.12 (0.13)	69.0
3.4 (0.110)	12.4 (0.81)	1.58 (0.06)	<lod< td=""><td>3.46 (0.49)</td><td><lod< td=""><td>20.8</td></lod<></td></lod<>	3.46 (0.49)	<lod< td=""><td>20.8</td></lod<>	20.8
	4250 (274) 14.4 (2.44) 3.4 (0.110)	IXN IXN-C 4250 151 (274) (1.86) 14.4 28.1 (2.44) (0.52) 3.4 12.4 (0.110) (0.81)	IXN IXN-C 8-PN 4250 151 4.46 (274) (1.86) (0.16) 14.4 28.1 5.93 (2.44) (0.52) (0.52) 3.4 12.4 1.58 (0.110) (0.81) (0.06)	IXN IXN-C 8-PN 6-PN 4250 151 4.46 75.5 (274) (1.86) (0.16) (11.4) 14.4 28.1 5.93 <lod< td=""> (2.44) (0.52) (0.52) 3.4 12.4 1.58 <lod< td=""> (0.110) (0.81) (0.06)</lod<></lod<>	IXN IXN-C 8-PN 6-PN XN 4250 151 4.46 75.5 6640 (274) (1.86) (0.16) (11.4) (78.2) 14.4 28.1 5.93 <lod< td=""> 17.9 (2.44) (0.52) (0.52) (0.49) 3.4 12.4 1.58 <lod< td=""> 3.46 (0.110) (0.81) (0.06) (0.49)</lod<></lod<>	IXN IXN-C 8-PN 6-PN XN XN-C 4250 151 4.46 75.5 6640 76.9 (274) (1.86) (0.16) (11.4) (78.2) (1.54) 14.4 28.1 5.93 <lod< td=""> 17.9 3.12 (2.44) (0.52) (0.52) (0.49) (0.13) 3.4 12.4 1.58 <lod< td=""> 3.46 <lod< td=""> (0.110) (0.81) (0.06) (0.49)</lod<></lod<></lod<>

guidelines for analytical methods of trace analysis (31). In future applications this method is applicable to measuring food samples, for freshness, food fraud and due to SIDA could be applied to biological samples as well. The validation results support that deuteration of PF is an easy pathway to producing labeled PFs and apply them to SIDA in quantifying PF in beer, hop tea and hops. The ongoing research regarding the bioactivities of these compounds provides a need for further analytical methods to be developed. The next step would be applying these methods to future clinical trial samples, such as blood, feces, and urine after consumption of beer or supplements. Reports of hop PF are already confirmed by urine analysis for

biomarkers or in animal studies although SIDA is yet to be utilized (12, 48).

DATA AVAILABILITY STATEMENT

The raw data supporting the conclusions of this article will be made available by the authors, without undue reservation.

AUTHOR CONTRIBUTIONS

LB and CU wrote the manuscript and did the MS and NMR analysis of compounds synthesized. LB and MR

designed the experiments for the microwave synthesis of the labeled XN, 6-PN, 8-PN, IXNC and XNC. LB carried out the synthetic work for all analytes excluding the stryker cat reaction of XN. LB and SS carried out the SIDA analysis including interpretation of the data. LB, SS, and CU carried out the preparative HPLC. SS performed the HPLC calibration graphs of XN-C and 8-PN. HR and CU designed the stryker cat synthesis of the xanthohumol and carried out the synthetic work for the DXN. All authors contributed to reading the article and proofing it before publication.

FUNDING

The European Foundation for Alcohol Research (ERAB) reference EA 18 69. Wissenschaftliche Station für Brauerei in München e.V.

REFERENCES

- Nikolic D, Van Breemen RB. Analytical methods for quantitation of prenylated flavonoids from hops. *Curr. Anal. Chem.* (2013) 9:71–85. doi: 10.2174/157341113804486554
- Venturelli S, Burkard M, Biendl M, Lauer UM, Frank J, Busch C. Prenylated chalcones and flavonoids for the prevention and treatment of cancer. *Nutrition*. (2016) 32:1171–8. doi: 10.1016/j.nut.2016.03.020
- Barrea L, Annunziata G, Muscogiuri G, Arnone A, Tenore GC, Colao A, et al. Could hop-derived bitter compounds improve glucose homeostasis by stimulating the secretion of GLP-1? *Crit Rev Food Sci Nutr.* (2019) 59:528– 35. doi: 10.1080/10408398.2017.1378168
- Bolton JL, Dunlap TL, Hajirahimkhan A, Mbachu O, Chen SN, Chadwick L, et al. The multiple biological targets of hops and bioactive compounds. *Chem. Res. Toxicol.* (2019) 32:222–33. doi: 10.1021/acs.chemrestox.8b00345
- Osorio-Paz I, Brunauer R, Alavez S. Beer and its non-alcoholic compounds in health and disease. Crit Rev Food Sci Nutr. (2020) 60:3492–505. doi: 10.1080/10408398.2019.1696278
- Milligan S, Kalita J, Pocock V, Heyerick A, Cooman L, Rong H, et al. Oestrogenic activity of the hop phyto-estrogen, 8-Prenylnaringenin. *Reproduction*. (2002) 123:235–42. doi: 10.1530/rep.0.1230235
- Venturelli S, Niessner H, Sinnberg T, Berger A, Burkard M, Urmann C, et al. 6and 8-prenylnaringenin, novel natural histone deacetylase inhibitors found in hops, exert antitumor activity on melanoma cells. *Cell Physiol Biochem.* (2018) 51:543–56. doi: 10.1159/000495275
- Chia-Leeson O, Riepl H, Urmann C. Biotransformation of the promising neuro-regenerative hop chalcone xanthohumol c. *Planta Med.* (2019) 85:277. doi: 10.1055/s-0039-3399972
- Kirchinger M, Bieler L, Tevini J, Vogl M, Haschke-Becher E, Felder TK, et al. Development and characterization of the neuroregenerative xanthohumol c/Hydroxypropyl-β-cyclodextrin complex suitable for parenteral administration. *Planta Med.* (2019) 85:1233–41. doi: 10.1055/a-1013-1276
- Oberbauer E, Urmann C, Steffenhagen C, Bieler L, Brunner D, Furtner T, et al. Chroman-like cyclic prenylflavonoids promote neuronal differentiation and neurite outgrowth and are neuroprotective. J Nutr Biochem. (2013) 24:1953–62. doi: 10.1016/j.jnutbio.2013.06.005
- Pichler C, Ferk F, Al-Serori H, Huber W, Jager W, Waldherr M, et al. Xanthohumol prevents DNA damage by dietary carcinogens: results of a human intervention trial. *Cancer Prev Res.* (2017) 10:153– 60. doi: 10.1158/1940-6207.Capr-15-0378
- 12. Paraiso IL, Revel JS, Choi J, Miranda CL, Lak P, Kioussi C, et al. Front cover: targeting the liver-Brain axis with hop-Derived flavonoids improves lipid

ACKNOWLEDGMENTS

Many thanks to Friedrich Ampenberger (Dipl.-Brmst.) for producing the hop free beer at TUM Research Center Weihenstephan for Brewing and Food Quality. Many thanks to Dr. Richard Hammerl and Dr. Oliver Frank for carrying out the NMR measurements at the chair of Food Chemistry and Molecular Sensory Science at the Technical University of Munich (TUM). Moreover, we thank the Barth Haas group for donating the CO_2 spent hop extract for xanthohumol purification. We would also like to thank Dr. Biendl from the Hallertauer Hopfenveredelungsgesellschaft m.B.H. to provide the hop product for the isolation of XN used in the synthesis of DXN.

SUPPLEMENTARY MATERIAL

The Supplementary Material for this article can be found online at: https://www.frontiersin.org/articles/10.3389/fnut.2020. 619921/full#supplementary-material

metabolism and cognitive performance in mice. *Mol Nutr Food Res.* (2020) 64:2070034. doi: 10.1002/mnfr.202070034

- Zhang Y, Bobe G, Revel JS, Rodrigues RR, Sharpton TJ, Fantacone ML, et al. Improvements in metabolic syndrome by xanthohumol derivatives are linked to altered gut microbiota and bile acid metabolism. *Mol Nutr Food Res.* (2020) 64:1900789. doi: 10.1002/mnfr.201900789
- 14. De Keukeleire J, Ooms G, Heyerick A, Roldan-Ruiz I, Van Bockstaele E, De Keukeleire. D. Formation and accumulation of α -Acids, β -Acids, desmethylxanthohumol, and xanthohumol during flowering of hops (Humulus lupulus l.). J Agric Food Chem. (2003) 51:4436-41. doi: 10.1021/jf034263z
- Kac J, Mlinari,č A, Umek A. HPTLC determination of xanthohumol in hops (humulus lupulus l.) and hop products. *J Planar Chromat.* (2006) 19:58– 61. doi: 10.1556/JPC.19.2006.1.10
- Moreira MM, Carvalho AM, Valente IM, Gonçalves LM, Rodrigues JA, Barros AA, et al. Novel application of square-Wave adsorptivestripping voltammetry for the determination of xanthohumol in spent hops. J Agric Food Chem. (2011) 59:7654–8. doi: 10.1021/jf20 11485
- Wyns C, Derycke L, Soenen B, Bolca S, Deforce D, Bracke M, et al. Production of monoclonal antibodies against hop-derived (Humulus lupulus l.) prenylflavonoids and the development of immunoassays. *Talanta*. (2011) 85:197–205. doi: 10.1016/j.talanta.2011.03.047
- Rychlik M, Asam S. Stabilisotopenverdünnungsanalysen zur quantifizierung organischer spurenkomponenten in der lebensmittelanalytik. Umweltwissenschaften und Schadstoff-Forschung. (2009) 21:470– 82. doi: 10.1007/s12302-009-0082-0
- Berwanger S, Zapp J, Becker H. 13 C-labelling of xanthohumol in hop cones (Humulus lupulus). Planta Med. (2005) 71:530–4. doi: 10.1055/s-2005-864154
- Berwanger S, Frank N, Knauft J, Becker H. Biosynthetic 14C-labelling of xanthohumol in hop cones. *Mol Nutr Food Res.* (2005) 49:857– 60. doi: 10.1002/mnfr.200500040
- Ellinwood DC, El-Mansy MF, Plagmann LS, Stevens JF, Maier CS, Gombart AF, et al. Total synthesis of [¹³C]₂₋, [¹³C]₃₋, and [¹³C]₅-isotopomers of xanthohumol, the principal prenylflavonoid from hops. *J Label Compd Rad.* (2017) 60:639–48. doi: 10.1002/jlcr.3571
- 22. Ellinwood D, Blakemore P. Total synthesis of a series of 13C-labeled isotopomers of xanthohumol: the primary prenylated chalcone in the inflorescence of hops plants. *Abstr Pap Am Chem Soc.* (2017) 253:1. doi: 10.1021/np070158y
- 23. Brestensky DM, Huseland DE, McGettigan C, Stryker JM. Simplified, "one-pot" procedure for the synthesis of [(Ph3P)CUH]6, a stable copper

hydride for conjugate reductions. *Tetrahedron Lett.* (1988) 29:3749–52. doi: 10.1016/S0040-4039(00)82105-X

- Lee DW, Yun J. Direct synthesis of stryker's reagent from a cu(II) salt. Tetrahedron Lett. (2005) 46:2037–9. doi: 10.1016/j.tetlet.2005.01.127
- Vaidyanathan S, Surber BW. Microwave mediated hydrogen deuterium exchange: a rapid synthesis of 2H-substituted benzimidazole. *Tetrahedron Lett.* (2005) 46:5195–7. doi: 10.1016/j.tetlet.2005.05.110
- 26. Fang Z, Zhou GC, Zheng SL, He GL, Li JL, He L, et al. Lithium chloride-catalyzed selective demethylation of aryl methyl ethers under microwave irradiation. J Mol Catal A-Chem. (2007) 274:16–23. doi: 10.1016/j.molcata.2007.04.013
- Urmann C, Riepl H. Semi-Synthetic approach leading to 8-Prenylnaringenin and 6-Prenylnaringenin: optimization of the microwave-assisted demethylation of xanthohumol using design of experiments. *Molecules*. (2020) 25:4007. doi: 10.3390/molecules25174007
- Stevens JF, Taylor AW, Deinzer ML. Quantitative analysis of xanthohumol and related prenylflavonoids in hops and beer by liquid chromatography-tandem mass spectrometry. J Chromatogr A. (1999) 832:97–107. doi: 10.1016/S0021-9673(98)01001-2
- Vogelgesang J, Hadrich J. Limits of detection, identification and determination: a statistical approach for practitioners. *Accredit Qual Assur.* (1998) 3:242–55. doi: 10.1007/s007690050234
- Mandel J. The Statistical Analysis of Experimental Data. New York, NY: Dover Publications (1964).
- Verlag M. Chemische Analytik; Nachweis-, Erfassungs- Und Bestimmungsgrenze; Ermittlung Unter Wiederholbedingungen Begriffe, Verfahren, Auswertung, Berlin: Auswertung, Beuth Verlag (2008).
- 32. Miles CO, Main L. The kinetics and mechanism, and the equilibrium position as a function of pH, of the isomerisation of naringin and the 4'-rhamnoglucoside of 2',4,4',6'-tetrahydroxychalcone. J Chem Soc Perkin Trans. (1988) 2:195–8. doi: 10.1039/P29880000195
- Atzrodt J, Derdau V, Fey T, Zimmermann J. The renaissance of H/D exchange. Angew. (2007) 46:7744–65. doi: 10.1002/anie.200700039
- Mahoney WS, Brestensky DM, Stryker JM. Selective hydride-mediated conjugate reduction of.alpha..beta.-unsaturated carbonyl compounds using [(Ph3P)CuH]6. J Am Chem Soc. (1988) 110:291–3. doi: 10.1021/ja00209a048
- Stevens JF, Ivancic M, Hsu VL, Deinzer ML. Prenylflavonoids from humulus lupulus. *Phytochemistry*. (1997) 44:1575– 85. doi: 10.1016/s0031-9422(96)00744-3
- 36. Wyns C, Bolca S, De Keukeleire D, Heyerick A. Development of a high-throughput IC/APCI-MS method for the determination of thirteen phytoestrogens including gut microbial metabolites in human urine and serum. J Chromatogr B Analyt Technol Biomed Life Sci. (2010) 878:949– 56. doi: 10.1016/j.jchromb.2010.02.022
- Sus N, Schlienz J, Calvo-Castro LA, Burkard M, Venturelli S, Busch C, et al. Validation of a rapid and sensitive reversed-phase liquid chromatographic method for the quantification of prenylated chalcones and flavanones in plasma and urine. NFSJ. (2018) 10:1–9. doi: 10.1016/j.nfs.2017.11.001
- Turowski M, Yamakawa N, Meller J, Kimata K, Ikegami T, Hosoya K, et al. Deuterium isotope effects on hydrophobic interactions: the importance of

dispersion interactions in the hydrophobic phase. J Am Chem Soc. (2003) 125:13836–49. doi: 10.1021/ja036006g

- Stevens JF, Page JE. Xanthohumol and related prenylflavonoids from hops and beer: to your good health! *Phytochemistry*. (2004) 65:1317– 30. doi: 10.1016/j.phytochem.2004.04.025
- Boronat A, Soldevila-Domenech N, Rodríguez-Morat, ó J, Martínez-Huélamo M, Lamuela-Raventós RM, de la Torre R. Beer phenolic composition of simple phenols, prenylated flavonoids and alkylresorcinols. *Molecules*. (2020) 25:2582. doi: 10.3390/molecules25112582
- Kamiński DM, Gaweda K, Arczewska M, Senczyna B, Gago, M. A kinetic study of xanthohumol cyclization to isoxanthohumol – a role of water. J Mol Struct. (2017) 1139:10–16. doi: 10.1016/j.molstruc.2017.03.027
- Magalhães PJ, Dostalek P, Cruz JM, Guido LF, Barros AA. The impact of a xanthohumol-Enriched hop product on the behavior of xanthohumol and isoxanthohumol in pale and dark beers: a Pilot scale approach. *J Inst Brew.* (2008) 114:246–56. doi: 10.1002/j.2050-0416.2008.tb00335.x
- 43. Stevens JF. Xanthohumol and Structurally Related Prenylflavonoids for Cancer Chemoprevention and Control. In: Pezzuto JM, Vang O, editors. Natural Products for Cancer Chemoprevention: Single Compounds and Combinations, Cham: Springer International Publishing (2020). p. 319–350.
- Pinto C, Cestero JJ, Rodríguez-Galdón B, Macías P. Xanthohumol, a prenylated flavonoid from hops (Humulus lupulus l.), protects rat tissues against oxidative damage after acute ethanol administration. *Toxicol.* (2014) 1:726–33. doi: 10.1016/j.toxrep.2014.09.004
- Intelmann D, Batram C, Kuhn C, Haseleu G, Meyerhof W, Hofmann T. Three tAS2R bitter taste receptors mediate the psychophysical responses to bitter compounds of hops (Humulus lupulus l.) and beer. *Chemosens Percept.* (2009) 2:118–32. doi: 10.1007/s12078-009-9049-1
- Schönberger C. *The Processing of Hops*. Sawston; Cambridge: Woodhead Publishing (2006). p.123–148. doi: 10.1533/9781845691738.123
- 47. Magalhães PJ, Guido LF, Cruz JM, Barros AA. Analysis of xanthohumol and isoxanthohumol in different hop products by liquid chromatography-diode array detection-electrospray ionization tandem mass spectrometry. *J Chrom A*. (2007) 1150:295–301. doi: 10.1016/j.chroma.2006.08.019
- Quifer-Rada P, Martínez-Huélamo M, Chiva-Blanch G, Jáuregui O, Estruch R, Lamuela-Raventós RM. Urinary isoxanthohumol is a specific and accurate biomarker of beer consumption. J Nutr. (2014) 144:484– 8. doi: 10.3945/jn.113.185199

Conflict of Interest: The authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

Copyright © 2020 Buckett, Schinko, Urmann, Riepl and Rychlik. This is an openaccess article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) and the copyright owner(s) are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.